

ヒトの唾を加えた溶液中のインフルエンザウイルスの SGNP による濃縮

| Q-PCR | | 唾液 (μl) | 生食 (μl) | SGNP 10 μl 添加 | Ct(CP) | Tm #1 | |
|--------|---------|----------------------|----------------------|--------------------------|--------|-------|-------|
| 15 HAU | samples | 100 | 400 | no SGNP | 22.94 | 83.14 | |
| | | | | SGNP | sup | 24.58 | 83.04 |
| | | | | | ppt | 18.82 | 83.20 |
| | | 0 | 500 | no SGNP | 21.55 | 82.99 | |
| | | | | SGNP | sup | 24.34 | 83.12 |
| | | | | | ppt | 18.86 | 83.09 |
| | control | p.c. | | no SGNP | 14.58 | 83.38 | |
| | | n.c. | | | 36.4 | 78.09 | |

